Surgery System

Surgic AP

OPERATION MANUAL

Rx Only MADE IN JAPAN

Thank you for purchasing the Surgic AP. This product is a unit used for a dental surgery, and other surgical operation, etc. Read this Operation Manual carefully before use for operation instructions, care and maintenance. Keep this Operation Manual for future reference.

■ Classifications of equipment

- Type of protection against electric shock :
 - Class I equipment
- Degree of protection against electric shock :
 - Type BF applied part 
- Method of sterilization or disinfection recommended by the manufacturer :
 - See 7. Sterilization
- Degree of protection against ingress of water as detailed in the current edition of IEC60529 :
 - Foot Control ... IPX8 (Protected against the effects of continuous immersion in water)
- Degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide :
 - Foot Control: Category AP Equipment 
- Mode of operation :
 - Intermittent operation

Indications for Use

The device is a control unit to drive contra angles and other attachments that are used with appropriate burs and drills for drilling, reaming, decorticating, and smoothing of bone and other bone related tissue in a variety of dental, oral surgical and other small surgical procedures.

Cautions for handling and operation

- Read these cautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
 WARNING	A hazard that could result in bodily injury or damage to the device if the safety instructions are not followed.
 CAUTION	A hazard that could result in light or moderate bodily injury or damage to the device if the safety instructions are not followed.
NOTICE	General information needed to operate the device safely.

WARNING


The system may present a possibility of malfunction when used in the presence of electromagnetic interference wave. Do not install the system in the vicinity of the device which emits magnetic waves. Turn off the power switch of the Control Unit of this system when an ultrasonic oscillation device or an electrode knife located in the vicinity is used.

CAUTION

- Surgic AP is intended for use in dental, oral surgical, and surgical procedures by qualified personnel.
- Patient safety is a priority.
- Read this Operation Manual before use, and fully understand the functions of each part for starting use.
- Inspect the operating status of the equipment before use, and use only after confirming that no abnormalities exist.
- Test run the product to ensure its correct operation prior to using it.
- If the product should ever malfunction (excessive vibration, noise, heat, etc) please turn it off immediately and return it to your Authorized Dealer for inspection.

⚠ CAUTION

- When the product is very frequently used please consider the maintenance of a small stock of replaceable parts.
- To avoid possible injury or product damage, ensure that the Motor has completely stopped before changing burs.
- Severe shock – Eg. Dropping the product – may cause damage.
- Do not operate the control unit without the irrigation pump.
- Do not use the Irrigation Pump if its package has been damaged or opened.
- Remove the irrigation pump after operation is stopped.
- Do not bend the Irrigation Tube while the water pump is operating. It could cause tube breakage.
- Never attempt to disassemble the Control Unit, the Foot Control or the Motor.
- Handpiece Attachments should be cleaned, lubricated and sterilized immediately after use.
- Do not lubricate the Motor. Oil could generate excessive heat and cause damage.
- The Control Unit and the Foot Control cannot be sterilized by any method.
- The Control Unit may be cleaned with a moist cloth. Disconnect the power supply before cleaning.
- Do not clean the Control Unit with any solvent solutions.
- Do not disconnect the motor cord from the motor.
- Be sure to dispose of the irrigation tube as medical waste after use.
- Grounding reliability can only be achieved when the equipment is connected to an equipment receptacle marked "Hospital Only" or "Hospital Grade".
- To prevent the possibility of getting an electric shock, connect to an electrical outlet with protective earth.
- The system functions normally in the environment where the temperature is at 0-40°C (32-104°F), humidity at 10-85% RH, atmospheric pressure at 700- 1060hPa, and no moisture condensation in the Control Unit. Use at outside of these limits may cause malfunction.
- The maximum surface temperature of the motor might reach 45°C by the situation used at atmosphere 25°C.
- Follow the table to the right for the motion time of motor and handpiece. Continuous use over a long time may cause the overheating of the handpiece, which in turn may result in an accident.
- Applied parts for patient and/or operator are/is motor and Handpiece.
- Surgic AP needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information.
- Portable and mobile RF communications equipment can affect Surgic AP. Do not use RF equipment in close proximity to the product.
- The use of ACCESSORIES, transducers and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of Surgic AP as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of Surgic AP.
- Surgic AP should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, Surgic AP should be observed to verify normal operation in the configuration in which it will be used.
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.


Continuous time	
ON	OFF
3 min.	10 min.

NOTICE

- To separate from Supply mains, Turn off the Power Switch after each use.
- See the Operation Manual for Handpiece on how to use it.
- For service requirements and spare parts please contact your dealer.
- The use of NSK genuine pre-sterilized, disposable Irrigation Tube Kit is recommended.
- Store and transportation the system in the place where the temperature is at -10-60°C (14-140°F), humidity at 10-85%RH, atmospheric pressure at 500-1060 hPa, and the system is not subject to air with dust, sulfur, or salinity.
- E Type Non-OPTIC handpiece can use with this product.


1.Package Contents

No.	Description	Quantity
1	Control Unit	1
2	Foot Control (With Cord) (2m)	1
3	AC Electrical Cord (2m)	1
4	Non-Optic Handpiece	1
5	Non-Optic Motor (with Motor Cord) (2m)	1
6	Irrigation Tube	5
7	Coolant Solution Hanger Post	1
8	Motor Cradle	1
9	Tube Holder	7
10	Internal Irrigation Nozzle	1

No.	Description	Quantity
11	Irrigation Pump	1
12	Y-Connector	1
13	E-type Spray Nozzle	1
14	Nozzle Cleaning Fine Wire	1
15	Nozzle Holder	1
16	Protection Plug	1
17	O-ring (Protection Plug)	2

2.Control Unit

Description of Operation

◆ Keys on the Unit

(1) Program key

This key is used to select any one of 8 available programs. Press [+] to ascend program numbers and [-] to descend the numbers. By pressing either [+] or [-] the numbers rotate continuously through all available programs.

(2) Speed key

This key is used to set the motor speed. Press [+] to increase speed by one step, and [-] to decrease speed by one step. When [+] is pressed and the speed setting reaches maximum or [-] is pressed and the speed setting reaches minimum, an audible intermittent beep sounds, and the speed cannot be changed any further.

(3) Torque key

This key is used to set the torque. Press [+] to increase torque by one step and [-] to decrease torque by one step. When [+] is pressed and the torque setting reaches maximum, or [-] is pressed and the torque setting reaches minimum, an audible intermittent beep sounds, and the torque cannot be changed any further. The range of torque setting steps vary according to the gear ratio selected to match the handpiece attachment in use.

(4) System key

This key is used to calibrate the handpiece attachment before use. To activate the automatic calibration mode.

(5) Gear Ratio key

This key is to select match ratio of the handpiece attachment, before use, to the unit. Press this key until the LCD display exhibits the correct gear ratio of the handpiece attachment.

(6) Coolant Flow key

This key is used to select the coolant solution flow volume. 5 flow volume rates are available for selection, plus the flow can be turned off.


(7) Forward/Reverse key

This key is used to change the rotational direction of the motor. Press this key once to change the rotational direction.

(8) Memory key

This key is used to memorize the program parameters set by the operator. Press this key for approx. 1 second to memorizes parameter. An audible beep confirms that new program parameters have been memorized.

◆ LCD display on the Control Unit

(1) Coolant Flow Level

Displays the selected coolant solution flow volume level. The selected Flow level is indicated by 1 to 4 levels of illuminated indicators. No light indicates the coolant solution flow is off.

(2) Program Number

Displays the selected program number.

(3) Gear Ratio

Displays the gear ratio of the handpiece.

(4) Forward/Reverse Indicator

Displays the rotational direction of the motor.

(5) Speed

Displays the speed of the motor.

(6) Torque

Displays the Torque of the motor.

(7) Torque Bar Graph

During operation displays an approximate percentage indication of the actual operating torque relevant to the preset maximum torque. When all bars illuminate, the operating torque is at maximum. When bars are half illuminated then the operating torque is approximately 50% of the preset torque.

CAUTION

The LCD display panel is produced from liquid crystal and should always be treated with care.

3. Foot Control

Speed Control Pedal

This pedal is used to start and stop the motor and to vary the speed during operation.

4. Installation

4-1 Connecting the Motor Cord

Face the ▲ mark on the Motor Cord plug upward then insert the plug into the Motor Cord jack on the Control Unit (Fig. 1). A click is heard when the motor cord plug is correctly inserted into the control unit. To disconnect the plug, pull back the lock joint, then disconnect the cord (Fig. 2).

Fig.1


Fig.2

4-2 Connecting the Handpiece

- Insert the handpiece until it reaches the end of Motor Insert and screw it in with pushing until it locks with the click sound. (Fig. 3)
- To disconnect, hold the motor and the handpiece, and pull out straightly.


Fig.3

4-3 Connecting the Foot Control

Face the screw on the Foot Control Plug downward then insert the plug into the Foot Control Cord Jack on the Control Unit. Secure the plug by fastening the lock nut. See Figs. 4 & 5.



Fig.4

Screw

Fig.5

4-4 Connecting the AC Electrical Power Cord

Align correctly then insert the electrical power cord into the power cord connection jack at the back of the Control Unit (Fig. 6).


Fig.6

4-5 Connecting the Irrigation Pump and the Irrigation Tube

- 1) Engage the hole of Irrigation Pump with the Drive Shaft of Control Unit. (Fig. 7)
- 2) Set the pawl of Irrigation Pump to the mounting hole and turn it clockwise until it stops.
- 3) Insert the irrigation tube into the groove on the reverse side with slack between pump and groove. (Fig. 8)

Fig.7

- To disconnect, turn the irrigation pump counterclockwise until it clicks.

CAUTION

- Before replacing the irrigation pump, wipe off excess water on pump and drive shaft. The wet drive shaft and rollers can be slippery and cause improper operation.
- When mounting the Irrigation Pump, wipe off dirt and debris from the rotation axis of the Control Unit with alcohol-immersed cotton swab or cloth.
- Insert the replacement irrigation pump into the drive shaft straight (slow, soft) to prevent damaging rollers in pump.
- Run the replaced irrigation pump about 10 seconds on largest setting of water volume before operation to adopt irrigation tube to new pump.

Fig.8

- 4) Attach the Irrigation Tube Connector on the side of Irrigation Pump connected to the groove to the other Irrigation Tube Connector and screw it in until it stops.


Fig.9

4-6 Mounting the Coolant Solution Bottle

Insert the Coolant Solution Bottle Hanger Post into the holder on the Control Unit. Place the bottle as shown in Fig. 10.


Fig.10

4-7 Insertion of the Irrigation Tube

- 1) Place the coolant solution bottle on the hanger post and insert the irrigation Tube Needle into the bottle cap. (Fig. 11)


Fig.11

- 2) Open the tube cap to supply air into the bottle. (Fig. 12)


Fig.12

CAUTION

Do not operate the irrigation pump if the tube is bent or the tube clamp is in the closed position. This could cause the tube to burst or slip out of the bottle.

4-8 Mounting the Irrigation Nozzle

3 irrigation methods are available; Internal, External or both of them, depending on tool and operation procedure.

- (1) Internal Irrigation Nozzle

Connect the Irrigation Tube to the Internal Irrigation Nozzle firmly. (Fig. 13)


Fig.13

- (2) External Irrigation Nozzle
Connect the Irrigation Tube to the External Irrigation Nozzle firmly. (Fig. 14)


Fig.14

- (3) Internal and External Irrigation together
 1) Connect the ends of the Y-Connector to the Internal Irrigation Nozzle and External Irrigation Nozzle respectively as detailed in procedures (1) and (2).
 2) Connect the Irrigation Tube to the Y-Connector firmly.


Fig.15

4-9 Attaching the Tube Holder

Use the tube holder to combine together the Motor Cord and the irrigation Tube. It is easier to insert Motor Cord first and, next, the irrigation Tube. (Fig. 16)


Fig.16

5. Operation

5-1 Programming the Motor Operation

The Control Unit can memorize 8 sets of programs. Each program includes the following functions which will be automatically performed when the appropriate program number is selected.

- Gear ratio of contra angle handpieces
- Direction of rotation
- Coolant solution flow
- Speed
- Torque upper limit

- (1) Turn on the power by pushing the Main Power Switch toward [-]. Whenever the main power switch is turned ON, program number 1 is always displayed.

- (2) Select a program number.

Press the [Program]  key on the control panel until the program number you require is displayed.

- (3) Selecting the gear ratio of the handpiece relevant to the program.

Press the [Gear Ratio]  key the gear ratio of the handpiece to be used is displayed.

- (4) Setting the speed. Set the speed by pressing the [Speed]  key.

- Each time this key is pressed the display changes to the next speed level. By pressing this key for more than 1 second brings the speed quickly to the next level until the speed display reaches its upper or lower limit.
- When the speed setting reaches the upper or the lower limit, an audible beep is heard and the speed setting cannot be changed any further.

- (5) Setting the torque upper limit. Set the torque upper limit by pressing the [Torque]  key on the unit control panel.

Power Switch Symbol Mark		
Function	OFF	ON

- Each time this key is pressed the display changes to the next torque level. By pressing this key for more than 1 second brings the torque quickly to the next level until the torque display reaches its upper or lower limit.
 - When the torque setting reaches the upper or the lower limit, an audible beep is heard and the torque cannot be changed any further.
- (6) Select the rate of coolant solution flow volume. Select the rate of the coolant solution flow volume by pressing the [Coolant Flow]  key.
- The rate of coolant solution flow volume has 5 flow rate steps plus "no coolant flow".
- (7) Memorize settings.
- After completing steps 1-6 press the [Memory]  key for more than 1 second until a long audible beep is heard. The long beep confirms that the programming is completed. If you hear a short audible beep when the [Memory]  key is first pressed please ignore this signal and keep the [Memory]  key depressed until a long beep is heard.

Repeat the above steps 1-7 to program any one of the 8 available programs.

5-2 Calibration Function

The resistance of a handpiece attachment against the rotation of the motor varies slightly depending on the handpiece model, its age and condition, the degree of wear on the handpiece gears, and so on. The Surgic AP unit incorporates an automatic function to recognize the level of the resistance of any handpiece attached to the motor, and to calibrate the motor to rotate the handpiece attachment to the specific speed and settings required.

- 1) Press the [Gear Ratio]  key and select the gear ratio of the attached handpiece.
- 2) Press the [System]  key for approx. 2 seconds until a long beep is heard. "CAL" is displayed.

CAUTION

This equipment is optimized to obtain the highest accuracy at a gear ratio of 1/20. When using another gear ratio, please note that the accuracy decreases with an increase in the ratio relative to 1/20.

- 3) Attach the handpiece to the motor and remove the bur. Press the [System]  key again. After a moment the motor will automatically start to run. After a short series of resistance diagnosis is completed, the display returns to normal display and the motor will automatically stop. Calibration of the handpiece is now completed.

CAUTION

- Because the handpiece AUTOMATICALLY starts to run, any bur must be removed from the handpiece before beginning the calibration process. Leaving a bur in the handpiece may cause harm to the operator. Care should be exercised not to ever add any load to a handpiece during calibration, because an incorrect diagnosis would result in incorrect control.
- If "FAIL" is displayed on the liquid crystal panel, check the mounting of each part and operate again. If "FAIL" is still displayed, contact your dealer.

Calibration fail example

- If nonconformities such as extreme abrasion of the gear and mechanical loss exceed the control limit, correct data cannot be obtained during calibration and the calibration will fail. In this case, please request repair.
- If it has not failed, calibration fails in the case where the bar makes contact with something during calibration.

5-3 Standard Operation

- 1) Turn on the Main Power Switch. The Control Unit is ready to perform the program memorized in program number 1.
- 2) Program is selected pushing the Program Key.
- 3) Verify the details of the program on the display.

- 4) Operating the motor. Step on the speed control pedal and the motor will start to run. Speed increases as the pedal is depressed. When the pedal is fully depressed the speed reaches the maximum set value.
- 5) During operation of the motor, when the drilling load reaches the programmed torque upper limit the integrated torque limiter automatically activates to prevent torque application excessive to the set requirement. When the torque limiter activates, the motor stops over beeping after 5 seconds.
- 6) Stopping the motor Release the Foot Control Pedal, and the motor will automatically stop.
- 7) To reverse direction of the motor (and bur) simply depress the Control Unit Forward / Reverse Key. A warning beep can be heard when the direction is in reverse mode.

6. Care and Maintenance

6-1 Protection Circuit

An electronic circuit breaker automatically functions to protect the motor and the control unit if the motor is ever overloaded. Power supply to the motor will automatically be terminated and the Error code will be displayed on the control unit.

◆ Resetting the Protection Circuit

To reset the protection circuit, release and then depress the speed control pedal.

6-2 Error Code

If an operational problem occurs the display shows the error code to allow an immediate problem diagnosis.

Error Code Display	Error	Cause of Error	Remedy
E0	System Error	• Memory failure.	Please request repair.
E1	Excessive Current Detected	• Extended use under heavy load. • Short circuit in the motor cord.	Electrical contact may be insufficient. Securely re-connect the motor cord. When an error cannot be eliminated, request repair.
E2	Excessive Current Detected	Motor Cord failure.	
E3	Motor Sensor Error	• Motor sensor failure (Hall IC). • Motor cord failure. • Ingress of water into a motor.	
E4	Unit Interior Overheating Error	• Overheating by extended use under heavy load. • Operation of the unit under an extremely high temperature.	Allow it to cool down before use. In order that heat is sufficiently radiated, periphery of the main unit should be well-ventilated wherever possible. When an error cannot be eliminated, request repair.
E5	Braking Error	Failure in the start / stop switch circuit. (Abnormal voltage)	When rotation and stop are repeated in short frequencies, a circuit may be activated which limits acceleration at start. Wait a few seconds and then use. When an error cannot be eliminated, request repair.
E6	Motor Rotation Failure Error	When the motor stops for more than 5 seconds after reaching the torque upper limit.	This is not a failure. It stops for safety reasons. When the error is reset, this can be used as it is.
		• Handpiece attachment failure. • Motor failure.	Please request repair.

Error Code Display	Error	Cause of Error	Remedy
E 7	Pump excess voltage	The pump fails.	An error code is displayed even when the irrigation tube is normal, please request repair.
E 8	Pump excess current		
E 9	Foot control abnormality	Connection failure of the connector.	Check the connection state of the connector.
		Failure of the inside of the Foot Control.	The error code is displayed even when the connector is normally connected, request repair.

6-3 Replacement of the Fuse

If the control unit does not function, check the fuses. To access the fuse box simply squeeze the fuse box lock located on the side of the control unit (Fig. 17). If the lock is too tight use a pointed tool to squeeze the lock. Use only the specified fuse to prevent a fire.

Voltage	Order code	Rating
120V	1200131510	T3.15AL 250V
230V	1200116010	T1.6AL 250V


Fig.17

6-4 How to replace the O-rings at the motor insert

- 1) Remove the worn O-rings with a pointed tool like a needle.
- 2) Insert new O-rings the O-ring grooves.

* Optional O-ring: Order code D0312074080


Fig.18

6-5 Maintenance of the Control Unit and Foot Control

If blood or saline solution is stained on the control unit or foot control, remove the power code, wipe off the Control Unit or Foot Control with a damp cloth, and wipe off with the alcohol-absorbed cloth.

6-6 Maintenance of the Handpiece

- Clean and lubricate after each use and/or before autoclaving.
- Follow the Operation Manual for each handpieces.

6-7 Cleaning of Irrigation Pump

Clean the Irrigation Pump after each patient as follows;

- 1) Half fill a cup with clean water.
- 2) Remove the Irrigation Tube Needle from the bottle and soak the tip of the needle in water.
- 3) Start the pump and clean the needle for approximately 10 seconds at the maximum irrigation volume.
- 4) Take the needle out of water and remove the Irrigation Pump from the Control Unit.
- 5) Sterilize by autoclave after blowing the Irrigation Tube dry by air.


Fig.19

CAUTION

Be sure to clean the Irrigation Pump before autoclave sterilization. If saline solution remains inside the Irrigation Pump, the solution will solidify and might cause irrigation failure.

7. Sterilization

- Autoclave sterilization is recommended.
- Sterilization is required first time you use and after each patient as noted below.
- Can be autoclaved : Handpiece, Motor (Motor Cord), Motor Cradle, Internal Irrigation Nozzle, Tube Holder, Protection Plug, Irrigation Pump

CAUTION

Do not autoclave any parts (Control Unit, Foot Control, AC Electrical Cord, Irrigation Tube, Y-Connector, Fuse) other than those that can be subjected to autoclave sterilization.

NOTICE

The Irrigation Pump does not need sterilization before its first use because it is shipped sterile.

■ Autoclave Procedure

- 1) Remove blood and debris from the handpiece.
- 2) Clean inside the handpiece; by using the spray lubricant (refer to "6-6 Maintenance of the Handpiece").
 - Attach the Protection Plug to the motor.
- 3) Clean the Irrigation Pump (refer to "6-7 Cleaning of Irrigation Pump").
- 4) Place those in autoclave pouch (not included in the package) and seal it.
- 5) Autoclave up to max 135°C in accordance with the tables below.


Fig.20

<Autoclavable parts OTHER THAN the Irrigation Pump>

Temperature	Sterilization time
121°C	20 min.
132°C	15 min.

* Sterilization at 121°C for more than 15 minutes is recommended by EN13060 or EN ISO17665-1.

<Irrigation Pump>

Type	Temperature	Sterilization time
Gravity Displacement	132°C	15 min.
Pre-Vacuum (Dynamic Air Removal)	132°C	4 min.
	134°C	3 min.

- 6) Keep the product in the autoclave pouch to keep it clean until you use it.

⚠ CAUTION

- Clean and lubricate the handpiece before autoclaving. Autoclaving a handpiece stained with blood or debris could cause damage to the handpiece.
- The Irrigation Pump can be autoclaved up to 20 times.
*Spare Irrigation Pump (Pack of 2) : Order Code Y1002196
- The cover for the Irrigation Tube Needle cannot be autoclaved. Remove the cover when sterilizing the Irrigation Pump (Fig. 21). In that case, be careful that the Irrigation Tube Needle does not pierce the autoclave pouch. A needle protruding out of the pouch may cause injury.
- Keep the product in suitable atmospheric pressure, temperature, humidity, ventilation, and sunlight. The air should be free from fine dust, salt and sulphur.
- Do not lubricate the motor.
- Do not disconnect the motor cord from the motor.
- Do not sterilize the product by autoclave sterilization at over 135°C. This causes product malfunction.
- If the sterilizer chamber temperature exceeds 135°C during the dry cycle then delete the dry cycle.
- The irrigation tube is a single use disposable type and cannot be autoclaved.
- You can hang Protection Plug on the Hanger Post/anywhere to prevent it from being lost. Do not hang the Protection Plug with motor. It may drop the motor and cause damage.


Fig.21

NOTICE

Protection Plug O-ring is consumable supply. If it is damaged, replace it with a new one.

Spare Protection Plug: Order Code 10001595

Protection Plug O-Ring: Order Code D0312457102

8. Specifications

< Control Unit >

Product Name	Surgic AP
Model	NE189
Power Supply Voltage	AC120V 50/60Hz AC230V 50/60Hz
Power Consumption	50VA
Max. Pump Output	75mL/min
Dimensions	W207 x D275 x H101mm
Weight	2.7kg

< Non Optic Motor >

Model	SG50MS
Speed Range	200~40,000 min ⁻¹
Input Voltage	DC30V
Dimensions	ø23 x L109.6mm (Without the motor cord)
Weight	287g (Including the motor cord)

9. Disposing Product

Please consult with dealer from whom you purchased it about waste disposal.

Symbols

TUV Rhineland of North America is a Nationally Recognized Testing Laboratory (NRTL) in the United States and is accredited by the Standards Council of Canada to certify electro-medical products with Canadian National Standards.

The product is designed not to become the ignition source in air and flammable anesthetic gas.
* Foot Control ONLY.

Protected against the effects of continuous immersion in dust and water.

This product is Autoclavable up to Max.135°C.

Dispose of this device and its accessories via methods approved for electronic device and in compliance with the Directive 2012/19/EU.

Type BF applied part.

Refer to the Operation Manual.

Marking on the outside of Equipment or Equipment parts that include RF transmitters or that apply RF electromagnetic energy for diagnosis or treatment.

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

Guidance and manufacturer's declaration - electromagnetic emissions		
Surgic AP is intended for use in the electromagnetic environment specified below. The customer or the user of the Surgic AP should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR11	Group 1	Surgic AP uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	class B	
Harmonic emissions IEC61000-3-2	class A	
Voltage fluctuations/ flicker emissions EN/IEC61000-3-3	Complies	Surgic AP is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration - electromagnetic immunity			
Surgic AP is intended for use in the electromagnetic environment specified below. The customer or the user of the Surgic AP should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC61000-4-2	± (2, 4, 6) kV contact ± (2, 4, 8) kV air	± (2, 4, 6) kV contact ± (2, 4, 8) kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC61000-4-4	± 2kV for power supply lines ± 1kV for input/output lines	± 2kV for power supply lines ± 1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	± 1kV differential mode ± 2kV common mode	± 1kV differential mode ± 2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5% Ut (>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 sec	<5% Ut (>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Surgic AP requires continued operation during power mains interruptions, it is recommended that the Surgic AP be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: Ut is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity			
Surgic AP is intended for use in the electromagnetic environment specified below. The customer or the user of the Surgic AP should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC61000-4-6	3Vrms 150 kHz to 80MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Surgic AP, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC61000-4-3	3V/m 80MHz to 2.5 GHz	3 V/m	Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P} \text{ 80MHz to 800MHz}$ $d = 2.3\sqrt{P} \text{ 800MHz to 2.5GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range ^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobiles radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Surgic AP is used exceeds the applicable RF compliance level above, the Surgic AP should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Surgic AP. b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3 V/m.			

Cables and accessories	Maximum length	Complies with
Motor with Motor Cord	2m	RF emissions, CISPR11
Foot Controller with Cord	2m	Harmonic emissions,
AC Electrical Cord	2m	Voltage fluctuations/ flicker emission,
		Electrostatic discharge(ESD),
		Surge
		Voltage dips, short interruptions and voltage variations on power supply input lines
		Power frequency (50/60Hz) magnetic field
		Conducted RF
		Radiated RF
		Class B/ Group 1
		IEC61000-3-2
		IEC61000-3-3
		IEC61000-4-2
		IEC61000-4-5
		IEC61000-4-11
		IEC61000-4-8
		IEC61000-4-6
		IEC61000-4-3

Recommended separation distances between portable and mobile RF communications equipment the Surgic AP.

The Surgic AP is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Surgic AP can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Surgic AP as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	80MHz to 800MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800MHz to 2.5GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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NSK America Corp. www.nskdental.com